

Laboratory Drying Oven 2m3

The **Alsto** 2M3 Drying Oven is ideal for high throughput applications, designed to dry mineral samples prior to sample preparation and analysis.

In addition to their ease of installation and maintenance, these digital-controlled ovens optimize laboratory floor space effectively. Minimizing the footprint and maximizing available lab space, when installed side by side.

Features include:

- Maximum 120 degrees C
- Full stainless steel internal lining
- High temperature resistance insulation
- Two Trolleys with high temperature castors
- 0.37 kW recirculating fan
- 3 off 4.0 kW Heating Elements
- Digital Temperature Control - 3 phase
- Over temperature protection
- Control box is **central at the front** so multiple ovens can be placed side by side easily



The **Alsto** 2M3 Drying Oven comes with two drying trolleys Each trolley shelf can hold 10 x Alsto Small Drying Trays
or 4 x Alsto Large Drying Trays